# *sinbo.*

SBS 4447, SBS 4448 BODY FAT SCALE

INSTRUCTION MANUAL

EN

RU

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## SINBO SBS 4447 / SBS 4448 BODY FAT SCALE INSTRUCTION MANUAL

### \* Feature:

1. It can be used for weigh weight as normal weight scales
2. Fat measuring instrument, can weigh weight, at the same time, plan moisture rate measured fat rate, muscle a bone mass and basic metabolism etc. body composition.

### \* Specification:

**Weighing scope/dividing value:** 5 kg-180 kg / 0.1 kg

**Fat rate scope/dividing value:** 5-50% / 0.1%

**Moisture rate scope/dividing value:** 20-70% / 0.1%

**Applicable height range:** 75-225 cm

**Applicable scope of age:** 10-85 years old

**Capacity:** 5 kg-180kg

**Unit:** kg, lb. st.

It can store 12 personal data

**Overload Indication:** Err; Low power Indication: Lo

**Tempered safety glass platform:** 280\*280\*5mm

**Display:** LCD and backlight

**Size of display:** 75\*30mm

**Battery:** 2\*AAA

**Voltage display:** The battery voltage is less than 2.4 V, shows "Lo" 2 seconds after the auto power-off

**The auto power-off:** no operation about 10 seconds will automatically be turned off

**Function keys:** UP, SET, DOWN

### \* Instruction Manual

#### 1. Parameter Enactment

a. Press Key **SET** to turn on (**every time turn on scale to set it**), User serial number (n=0-11) flicker, press key **▲** or **▼** (increase by degrees or degression Number cycle) to choose the number which you want (this item can supply 12, using 0-11 expressing), then press **SET** to confirm. If there is no operation within 3 seconds, system will confirm automatically current flicker number is the current user measure data.

b. After user number enactment, gender sign start flicker. Now press **▲** or **▼** to choose gender, Press **SET** to confirm, if there is no operation, it will be turn off automatically within 20 seconds.

c. After gender enactment, age will start flicker. Now press **▲** or **▼** to adjust age, press

**SET** to conform. If there is no operation, it will turn off automatically within 20 seconds.

**d.** After age enactment, height will start flicker. Now press ▲ or ▼ to adjust height, press **SET** to conform. If there is no operation, it will turn off automatically within 20 seconds.

**e.** After All the data enactment, showing "0.0",mean you can start to weigh fat

## **2. Measure**

### **1. Measure weight**

**a.** Under the turn off, please move on reposefully the scale, screen will show your weight, after 3 times, it will be locked.

**b.** If there is no operation, it will be turn off automatically within 10 second.

### **2. Measure body fat, water, muscle, bone, calorie, weight, BMI**

Please bare feet stand steadily on the scale ,feet with two electrode contact , body don't swing .The screen will show your weight value , stability weight data flashing three time then be locked , it will indicate fat, water, muscle, bone, calorie, weight , each measurement value will be displayed for two seconds , alternating cycle three times after the automatic shutdown .

### **Usage Notes:**

**1.** Ensure that the scale is placed on a hard & level surface; otherwise the result will be inaccurate if on the carpet or uneven ground.

**2.** Keep scale surface dry to avoid accident by slipping.

**3.** When weighing please take off the socks, shoes and clothes, at relatively fixed time before meal.

**4.** In order to measure precision, not recommend to weight in 2 hours after getting up.

**5.** Please put the scale in a dry place to protect the internal electronic components.

**6.** Please use micro cloth to clean the scale, not immersed in water.

**7.** Please don't remove the scale except replace the battery.

**8.** Please don't use heavy percussion or falls from a height.

**9.** Please don't persistent put heavy on the scale to avoid damage to the high precision sensor.

### **CORRECT DISPOSAL OF THIS PRODUCT**

The symbol on the unit, packing materials and/or documentations means used electrical and electronic units and batteries should not be toss in the garbage with ordinary household garbage. These units should be pass to special receiving point.

For additional information about actual system of the garbage collection address to the local authority.

Valid utilization will help to save valuable resources and avoid negative work on the public health and environment which happens with incorrect using garbage.

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